QuickSee *Free Pro*
KERATOMETRY

QuickSee *Free*
QuickSee *Free Pro*
KERATOMETRY

English

Español

Français

Deutsch

Italiano

Nederlands

Dansk

Português

Svenska

Norsk

Türkçe

Lehçe

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User and safety instructions

N: TF-002-AN2-001 2025-12-17 Rev F

English

User and safety instructions

Congratulations on purchasing the PlenOptika QuickSee Free Pro Keratometry or the QuickSee Free! This booklet will help you make the most of your experience with your QuickSee Free medical device.

Intended use

QuickSee Free medical devices are intended to measure the refractive power of the eye by measuring light reflected from the eye using the optical techniques of wavefront aberrometry and keratometry (optional, only in QuickSee Free Pro Keratometry model).

Indications for use

The QuickSee Free medical devices are indicated for use by healthcare professionals or under the direction of a healthcare professional to screen for or measure an individual's refractive errors associated with poor vision.

Contraindications

The refractive error estimate is contraindicated in the following cases:

- Ocular drugs interfering with the visual system may affect the results.
- A history of eye surgery or the presence of mature cataracts may also affect the results.

Clinical benefits

1. Improved management of refractive errors (more than 2 billion people suffer from uncorrected refractive error according to WHO data). The portability, usability and features of the technology allow the eyecare professionals to be more efficient and enable minimally trained personnel to participate in the refraction process.
2. QuickSee Free medical devices support the intended use of the device and derived model, which is to produce an estimate of the refractive status of an eye.
3. QuickSee Free medical devices enable the intended use of the device and derived model, which is to provide estimates of sphere, cylinder, and axis to assist and aid the eye-care professional in evaluating refractive errors and prescribing vision correction for selected subjects. In addition, QuickSee Free Pro Keratometry also provides keratometry readings.

Incident reporting

Incidents with medical devices should be reported to the competent authority of the country in which the user and/or patient is established. If you have any doubt about the process, contact your local distributor or PlenOptika.

By default, PlenOptika will report incidents to The Spanish Agency of Medicines and Medical Devices (AEMPS). Incident reporting is an essential element through which AEMPS records reported incidents in order to identify any problems or defects related to medical devices. This registry enables the association of cases and the generation of indications of product failures, as well as the assessment of the impact of the events and, if necessary, the rapid adoption of the relevant corrective measures.

If you suspect an incident involving a medical device, you can report it using the electronic form found at the link below. To do so, choose one of the two available forms: healthcare professional or patient. More information on what, how, who, and why to report can be found in the link below:

<https://notificaps.aemps.es/enviotelematico/notificaps/notifica/inicio.do>

Device description

QuickSee Free Pro Keratometry is a hand-held wavefront aberrometer/keratometer with a monocular form factor. Operation involves asking the subject look at a distant object through the device and aligning their eye with the red dot and pupil camera displayed on the screen. The device provides a measurement of a subject's refractive error as well as keratometry values.


QuickSee Free is similar in all aspects of QuickSee Free Pro Keratometry but does not include the Keratometry functionality.

Box contents

The shipped package contains:

- QuickSee Free Pro Keratometry or QuickSee Free medical device
- Device case
- Medical grade USB power adapter and 4 international plug adapters
- Power cord
- Fogging lenses (2 pieces, detachable)
- Eyecup (detachable)
- Bluetooth printer
- Power cord for bluetooth printer
- QuickStart guide
- Occluder
- Wrist strap
- Supplemental materials

Terminology

Basic operation

Basic care

1. The QuickSee Free medical devices are precision instruments. **Do not** drop or subject it to physical shock. If the device is dropped or suffers a physical shock, functionality may be affected. To ensure proper functionality of the device after any non-negligible physical shock, the device needs to be inspected by PlenOptika.
2. The QuickSee Free medical devices are **not waterproof**. If you accidentally drop it into liquid or liquid is spilled on the device, contact PlenOptika Support immediately. If minor water droplets are found on the device, wipe off with a soft dry cloth. Please wipe off drops of liquid like saliva and water on LCD panel surface immediately.
3. **Do not** leave QuickSee Free in excessive heat, such as in a car or in direct sunlight. High temperatures can cause the device to malfunction.
4. QuickSee Free contains precision electronic circuitry. Never attempt to disassemble the device or remove the battery; this may result in bodily harm and will void any warranty and must be always performed by Service Personnel.
5. **Do not** use cleaners that contain organic solvents to clean the QuickSee Free's body, LCD screen, view ports, or back viewing lens. Using such solvents may damage the device.
6. Use a microfiber cloth for cleaning. For stubborn stains, wet the cloth or use a non-abrasive cleaner. Always keep the front and back view ports and fogging lens in the viewport free of fingerprints, dirt, and dust.
7. When plugging the power connector into QuickSee Free, **DO NOT** force the power connector into the device as this may cause damage and will not be covered under warranty.
8. **Do not** push or scratch the touchscreen display surface with anything stiff or rigid. Such an instrument may damage the device. A stylus is acceptable for use with the device if desired.
9. **Do not** store QuickSee Free in its case with the power cord plugged in. This could cause damage to both the power cord and the device.
10. The device is **not** intended to be used in conjunction with any flammable gases or liquids.
11. Electrically powered medical devices connected to the electrical network in which the QuickSee Free device is connected should meet EMC standards.
12. The device is **not** intended to be used in intense electromagnetic environments. Check EMC specifications of the device in this document.
13. The applied parts of the QuickSee Free medical devices (i.e., the parts of the medical devices that in normal use necessarily come into physical con-

tact with the patient to perform its function) are made of appropriate materials to minimize the probability of patient-to-patient contamination. It is recommended that these parts are regularly cleaned and disinfected as per normal procedure for optometric equipment.

14. The device contains an IEC 60825-1: 2014 certified Class 1 Laser Device (780 nm laser radiation) and 850 nm LED illumination ring sources emitting within the Acceptable Emission Limits during normal use and fault conditions.

Failure to observe these warnings and recommendations may void any warranty on the QuickSee Free medical devices.


Warnings

1. Changes or modifications not expressly approved by PlenOptika could void the warranty.
2. The supplied power cord (white color) is optimized to charge the QuickSee Free medical device. Replacement or substitution may only be done with a CSA or UL approved medical-grade power cord rated for at least 5V and 3.2 amp.
3. The supplied power supply/charger may not be substituted and shall only be replaced by an approved unit supplied by PlenOptika.
4. If the power cord or external power supply becomes damaged or its integrity is suspect, or the building power receptacle outlet is ungrounded, then QuickSee Free medical devices shall only be used with the internal battery power until replacement parts or a grounded power receptacle become available.
5. The power cord for the Bluetooth printer (black color) is not optimized to charge the QuickSee Free medical device, being exclusive for the Bluetooth printer. Likewise, the power cord for the QuickSee Free medical device (white color) shall be used exclusively for the medical device.
6. QuickSee Free medical devices contain a series of magnets in the device itself and in the detachable fogging lenses. The magnetic field emitted could cause minimal interference with susceptible implantable devices such as pacemakers, and it is recommended that QuickSee Free medical devices are not placed at a distance less than 5 cm from the patient's chest. Although operators should be aware of this recommendation, it should be noted that QuickSee Free medical devices are safe when used under normal operation conditions following the user and safety instructions.

Getting started

Charging


1. Locate the QuickSee Free charging port on the bottom side of the handle (see figure at right)
2. Connect the white power cord with supplied power adapter and plug into the charging port.
3. Plug the power adapter into an available AC wall outlet to charge the device; a purple light indicates the device is charging, and a blue light indicates the device is fully charged.
4. The unit can be used while charging. Ensure unit is charged fully before first use by charging for 8 hours initially.

Warning: Charging port specifications: 5V-3.2A DC Input. Use provided medical grade power cord.

Powering QuickSee Free on and off

1. To turn QuickSee Free medical device on, press the power button, located on the top side of the device (see figure below), for half a second and release. The device will power on and the display will turn on. Wait until the device finishes the startup procedure (around 45 seconds).
2. To turn the device off, access the pull-up menu (see how to access the pull-up menu on page 8) and select "Power off" by sliding the control from left to right.
3. Pressing and releasing the power button briefly will put the device into Pow-

er-Saving Mode. In Power-Saving Mode, the screen is turned off and less battery is consumed. The device is instantly operational when power button is pressed.


- If the QuickSee Free medical device display freezes and is unresponsive to touch, you will need to perform a hard system reset by holding the power button for 7 seconds until the screen goes black and then releasing the button.
- Low battery: When the battery level is very low, a warning message will appear. If charging is not initiated soon after, the device will automatically shut down.

User operation

Home screen


When the QuickSee Free medical device is powered on, the Main Menu or “Home Screen” will appear. If the device’s Bluetooth connection is enabled, the Bluetooth icon will be visible in the top right corner. Bluetooth is enabled by default when the device powers on and can be turned off in the General Settings. The icon will change if the device is in discoverable mode or connected, and will not appear if the Bluetooth connection is turned off.

- To begin a measurement, press OD or OS icon to start measuring the right or left eye, respectively.
- The drop-down or pull-up or menus (see page 8) can be accessed by sliding your finger from the top or bottom of the Home Screen, respectively.


Measuring a patient

- On Home Screen, select the OD or OS icon to start measuring the right or left eye, respectively.
- Bring the device to the eye you want to measure while looking at the live video of the patient’s eye. Ensure that the two contact points of the device are touching the patient’s face.

Patient and operator positions during measurement:


- The patient looks through the device at a distant target.
 - The operator stands at patient’s side and sees the patient’s eye on the display.
- Instruct the patient to hold the occluder over their other eye while keeping it open.
 - Align the red spot with the center of the pupil. The background of the display and the play icon will both turn green when the spot is properly aligned to the pupil.
 - Ensure that you are not in the patient’s line of sight, which will enable the patient to see the distant target and reduce patient’s accommodation. Note that the display can be rotated to give you a view of the display.
 - When the red spot is aligned with the patient’s pupil (the background turns green), press the play

icon on the screen or the start button in the handle to start a measurement. During the measurement, ask the patient to:

- a. focus on the distant target instead of the red dot, and
- b. blink at the start of the measurement and a few times during the measurement.



7. After the acquisition starts, a dotted line will appear around the pupil's image indicating the acquisition progress.
8. After 10 seconds (default value), the patient's estimated prescription for the selected eye will be displayed on the Results Screen.

Results screen


After measurements are taken, the Results Screen will show the estimated prescriptions. (NOTE: keratometry results will be shown only on QuickSee Free Pro Keratometry).

1. Swipe horizontally to navigate between eyes and vertically to navigate between measurements for the selected eye (if more than one measurement was performed).
2. To re-measure an eye: navigate to the specific eye, press the "+" button, and select OD or OS to start a new measurement of the same patient.


3. QuickSee Free Pro Keratometry only: Press the  symbol to view advanced refraction metrics. The advanced refraction metrics that can be displayed on QuickSee Free Pro Keratometry include:

- Pupil image
- Keratometry results
- Pupil size
- Spot diagram
- Zernike coefficients
- Wavefront map
- Simulated image, uncorrected
- Simulated image, corrected


Sample advanced refraction metrics screens

4. Measurements are saved for the current patient as soon as they are taken. Press the “+” button and slide the slider button from left to right to start a new patient record.



Drop-down menu

From any screen, slide the slider button from the top of the screen to the bottom. This will allow you to see the date, time, and battery level, to adjust the screen brightness, and to activate Bluetooth (Bluetooth icon) or the night-time mode (moon icon).

Pull-up menu

From any screen, slide the slider button from the bottom of the screen to the top. This will allow you to access settings, help, or about menus. Additionally, you can turn off the device from this screen.

Settings options

General

Date/Time
Language
Bluetooth
Printer
Screen timeout
Shutdown time

Prescription

Cylinder Notation
Diopters precision
Axis precision
Eye notation
K unit (Pro only)
K index (Pro only)

Measurement

Acquisition Time
Usage reminders at boot
Measurement tips
Measurement mode (Pro only)
Advanced refraction metrics (Pro only)

Reset

Erase patient data
Reset Bluetooth connection
Factory reset

Data management and software upgrades

Software upgrades and access to internal data stored in the QuickSee Free medical devices can be performed only with the QuickSee Free Companion App software provided by PlenOptika. All the information related to the QuickSee Free Companion App software (including download links) can be found at <http://plenoptika.com/user-manuals>

Minimum hardware requirements for the QuickSee Free Companion App

Connectivity (any operating system and hardware):

- Bluetooth 5.0 or higher (Android and Windows); Bluetooth 5.2 or higher (macOS).
- Internet functionalities require an active internet connection (cable or wireless) with a download speed of at least 3 Mbps.

Hardware requirements (any operating system):

- Storage: At least 500 MB of available disk space.
- At least 2 GB of RAM Memory.
- Any microprocessor capable of running one of the following operating systems.

Operating systems:

- Desktop computers running Windows 10 or 11, version 21H2 or later, Architecture x86_64. Note: PlenOptika will no longer provide support for Windows 10 once Microsoft no longer provides official support for this operating system.
- Desktop computers running macOS, Sonoma version 14.0 or later.
- Tablet or smartphone running Android 11 (Red Velvet Cake) or higher.

Data security

QuickSee Free medical devices and accessories have been designed taking into account cybersecurity risks, providing protection against unauthorized access, and are inherently safe. Follow these recommendations for enhanced protection:



- Keep the QuickSee Free medical device within your physical control.
- QuickSee Free medical devices have been designed to connect only to other devices using PlenOptika proprietary software (i.e., QuickSee Free Companion App). Do not try to connect QuickSee Free medical devices to other devices through unauthorized third-party software.
- Keep the QuickSee Free medical device software updated. Software updates may include security updates, such as patches for new cybersecurity risks.

Connecting to the QuickSee Free Companion App



There are two ways to connect the QuickSee Free medical device to the Companion App:

1. Connecting to the Companion App for the first time:

a. Put the QuickSee Free medical device in discoverable mode (Bluetooth). The device is in discoverable mode on initial startup by default, indicated by the Bluetooth searching  symbol. If the QuickSee Free medical device is not in discoverable mode, turn the Bluetooth off, and then on again, from the drop-down menu.

b. From the Companion App, turn on the Bluetooth of the PC, smartphone, or tablet (if it is off) and scan available devices. Then connect to the desired QuickSee Free medical device, identified by its serial number.


c. If it is the first time connecting to that QuickSee Free medical device from that device (PC, smartphone, or tablet), a random security code will be requested that has to be introduced in the Companion App within a minute. After submitting the correct security code, the Companion App will connect to the QuickSee Free medical device. If the security code was not requested, the Companion App will directly connect to the QuickSee Free medical device. For Windows systems, it is recommended to first pair the QuickSee Free medical device through Windows' Bluetooth settings.

2. Reconnecting the Companion App to a previously connected QuickSee Free medical device:
 - a. Turn the Bluetooth on from the drop-down menu of the QuickSee Free medical device or put the QuickSee Free medical device in discoverable mode (Bluetooth). Discoverable mode is indicated by the Bluetooth searching symbol. If the device is not in discoverable mode, turn the Bluetooth off, and then on again, from the drop-down menu.
 - b. From the Companion App, turn on the Bluetooth of the PC, smartphone, or tablet (if it is off) and connect directly to an already known QuickSee Free medical device identified by its serial number or by a custom alias previously set by the user.


When the QuickSee Free medical device is connected to the Companion App, the QuickSee Free medical device will show a link symbol and the Companion App will display the data from the device.

Connecting to the Bluetooth printer


Each QuickSee Free medical device is shipped with a pre-paired Bluetooth printer.

To print measurement results, turn on the paired Bluetooth printer with a long press on its power button. When results are available, you can print them from the Results Screen by pressing "+", and then pressing the printer symbol.

You can select the Bluetooth printer by accessing the settings, through the pull-up menu, selecting “General”, and then pressing “Printer”. The screen will subsequently show that nearby devices are being scanned, where you will have to select the “Bluetooth Printer” identified by its MAC address.

If you do not have a Bluetooth printer selected or the selected printer is not available, a warning message will appear.

Technical information

Technical description

The autorefractometer system is based on wavefront aberrometry, a technology currently used in benchtop wavefront aberrometers to measure low- and high-order refractive (wavefront) errors and to guide LASIK surgery. The specific implementation of this technique in the devices enables the user to obtain an open-view autorefractometer.

The keratometry system included in the QuickSee Free Pro Keratometry medical device determines corneal curvature by analyzing reflections of an illumination source off the cornea.

Measurement environment

The QuickSee Free medical devices can be used indoors or outdoors. Be sure to eliminate or block any sources of sunlight and/or light that may shine onto the patient's eyes during the measurement.

The device is designed to work in the following environmental conditions:

- Ambient temperature range: -10 – 40 °C (14 – 104 °F)
- Relative humidity range: 30 – 90%
- Atmospheric pressure: 800 – 1060hpa

NOTE: If the device is used continuously in environmental conditions in which the temperature is 40°C, some parts of the device may surpass 41°C. During such situations, the expected maximum temperature of the following parts (as determined by IEC 60601-1 testing) are:

- Eyecup: Maximum expected temperature < 43°C

Under these conditions, the device does not represent a safety risk for the patient, operator, or environment. However, it is recommended to reduce the temperature of the device in order to avoid an unpleasant experience for the patient. To accomplish this, we recommend:

- Turning off the device until the temperature of these parts is reduced.
- Avoiding exposure to a heat source (e.g., sunlight).

NOTE: The measurement time is 10 seconds (according to IEC 60601-1 the exposure limit is 60°C for 60 seconds or less).

NOTE: A minimum pupil size of 2 mm is recommended for measurements to be taken. It may be difficult to take a measurement in an environment with a bright light source that causes the pupil to constrict smaller than 2 mm. For best results, subject's pupils should be between 2 mm and 8 mm.

Transport and storage conditions

The QuickSee Free medical devices can be transported and stored in the following conditions:

- Ambient temperature range: -5 – 50 °C (23 – 122 °F)
- Relative humidity range: less than 80%
- Atmospheric pressure: 800 – 1060hpa

Environmental protection

QuickSee Free medical devices contain Lithium-ion cells and other electronic materials. Please contact your local government authorities and consult with relevant laws for proper electronic disposal methods or recycling at device end of life.

Servicing and technical documents

This documentation is part of the technical file of the QuickSee Free medical devices. PlenOptika will make available on request circuit diagrams, component part lists, descriptions, or other information that will assist PlenOptika authorized Service Personnel to repair the QuickSee Free medical devices.

EMC information

The QuickSee Free medical devices conform to the most widely accepted EMC standard, the UNE-EN 60601-1-2:2015. The device has been tested and found to comply with the limits for Group 1, Class A medical device (RF Emission CISPR 11), pursuant to Part 15 of the FCC Rules.

Electromagnetic compatibility

The QuickSee Free medical devices are suitable for use in all establishments other than domestic settings and may be used in domestic settings and those connected to the public low-voltage power supply network that supplies buildings used for domestic purposes, provided the following warnings are heeded.

- **Warning:** This equipment is intended for use by healthcare professionals or under the supervision of healthcare professionals only. This equipment may cause radio interference or may disrupt the operation of nearby equipment. It may be necessary to take mitigation measures, such as re-orienting or relocating the QuickSee Free or shielding the location.
- **Warning:** MEDICAL ELECTRICAL EQUIPMENT needs special precautions regarding EMC and needs to be installed and put into service according to the EMC information provided in the accompanying documents.
- **Warning:** Connectors should not be touched, and connections should not be made to these connectors unless ESD precautionary procedures are used.
- **Warning:** ESD protection mark  may be used in lieu of protective covers.
- **Warning:** Portable and mobile RF communications equipment can affect medical electrical equipment.
- **Warning:** The use of accessories and cables other than those supplied with the instrument, except for cables sold by the manufacturer of the equipment as replacement parts, may result in increased emission or decreased immunity of the equipment or system.
- **Warning:** The emissions characteristics of this equipment make it suitable for use in industrial areas and hospitals (CISPR 11 class A). If it is used in a residential environment (for which CISPR 11 class B is normally required) this equipment might not offer adequate protection to radio-frequency communication services.

FCC Compliance Statement

QuickSee Free medical devices comply with Part 15 of the FCC Rules. Operation is subject to the following two conditions: (1) this device may not cause harmful interference and (2) this device must accept any interference received, including interference that may cause undesired operation.

Technical specifications

Autorefractor

Pupil size	2 to 8 mm
Intended patient population	>= 3 years
Accommodation control	Open view, Fogging lens (optional)
Cycloplegia requirement	None
Dilation Requirement	None
Illumination requirements	None
Acquisition time	5, 10 seconds
Spherical Range	-13D to +10D, increments of 0.01D, 0.125D, 0.25D
Cylindrical Range	-8D to +8D, increments of 0.01D, 0.125D, 0.25D
Axial range	0-180°, increments of 1°, 5°, 10°
Base Technology	Wavefront aberrometry

Keratometer (QuickSee Free Pro Keratometry only)

Radius of curvature	5 to 12 mm in 0.01 mm increments
Corneal astigmatism	0 D to +/- 8 D (R 5 to 12mm)
Center	Diameter 3.2 mm at Radius of curvature 8 mm
Pupil size measurement	Displayed in advanced refraction metrics, Companion App, and printed results

Bluetooth and Wi-Fi specifications

Frequency Band BT	2.4 to 2.485 GHz
Frequency Band Wi-Fi	2.4 GHz and 5 GHz
Wireless Radio Model	Compliant with FCC
Antenna Type	Chip antenna, FHSS Modulation
Dataflow	Bidirectional

General specifications

Measurement modes (QuickSee Free Pro Keratometry only)	Simultaneous AR and K, AR only, and K only
Display properties	2.4-inch LCD, capacitive touch screen, readable outdoors, true color (65,536 colors)
Charger properties	USB-C medical grade wall adapter, AC 100 to 240V, 50/60 Hz
Battery	6 hours continuous use +/- 1 hour (10,000 mAh Li-ion); 3 hours charge time (5% - 75%); 5 hours charge time (0% - 100%); IEC 62133-2:2017 certified
Calibration	Factory calibrated; no field calibration needed
Measurement capacity	Measurement storage capacity: > 10,000 measurements
Regulatory classification (medical device)	Class I FDA (USA), Class IIa Product CE, MDR Compliant, Class IIa, UKCA, MDR 2002 Compliant
Laser safety	Class 1, IEC 60825-1:2014 certified
Weight	< 750 grams / 1.65 lbs
Dimensions	30 cm (H) x 5.5 cm (W), 18 cm (L)

Safety symbols

Symbol	Description
	Direct Current (DC) Power Supply.
	Product contains primary and secondary Lithium Ion batteries. Definition: Do not dispose the device and its batteries in standard trash containers. Proper disposal of Device and its components is required. Consult your nearest recycling center for proper disposal instructions.
	Manufacturer address.
	Manufacturing date.
	Consult Instructions for Use / User Manual.
	Type B Applied Part.
	Standby and on / off functionalities: <ul style="list-style-type: none"> • If the device is off, a short press on this button will turn on the device. • If the device is on, a short press will put the device in Standby mode. • If the device is in Standby mode, a short press will put the device in Standard mode. • If the device is on or in Standby mode, a long press (> 7 seconds) will force shutdown.
	Serial number of the device.
	Packaging, storage, and shipment. Indicates the correct upright position of the package.
	Packaging, storage, and shipment. The package has to be protected from moisture.
	Packaging, storage, and shipment. The package has to be handled carefully to avoid any damage.
	Packaging, storage, and shipment. Indicates the range of atmospheric pressure to which the package can be safely exposed.
	Packaging, storage, and shipment. Indicates the range of humidity to which the package can be safely exposed.
	Packaging, storage, and shipment. Indicates the range of temperature to which the package can be safely exposed.
	FCC marking. Indicates that the electronic device is certified and the electromagnetic interference from the device is under the limits that are approved by Federal Communications Commission.
	CE Marking Certification. Indicates that the product has met European Union consumer safety, health, and environmental requirements.
	UK Conformity Assessed marking is a conformity mark that indicates conformity with the applicable requirements for products sold within Great Britain.
	UDI (Unique Device Identifier) of the device.
	Medical Device.

Contact PlenOptika

If you have any questions about QuickSee Free medical devices or any other PlenOptika product, please contact customer support.

Email: support@plenoptika.com

Address:  PlenOptika Europe, S.L.U.
Parque Científico de Madrid
Calle Faraday 7, #2.06
Madrid, 28049 Spain

US office: PlenOptika, Inc.
955 Massachusetts Ave #339
Cambridge, MA 02139
USA

UKRP: Grafton Optical Company Limited
Unit 7 Riverpark Industrial Estate Billet Lane
Berkhamsted Hertfordshire
England, HP41HL
UK

Visit the PlenOptika website (www.plenoptika.com) for detailed guides, tips, and information on accessories and optional software features of QuickSee Free medical devices.

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